



Medisphere introduces a new line in Wound Care Support Surfaces - 1000

The reality of financial considerations is always a major factor in determining what therapies can be provided.

Medisphere is pleased to offer the Tradewind 1000 Mattress overlay system a versatile, affordable true alternating support surface as a solution that takes both the Patient and the bottom line into consideration .

Patients Assessed at Low to Medium Risk

www.medispherehealthcare.com

Technical Specifications

Tradewind 1000 Overlay System

Complete System Order Code: Med/1000/ovl

Pump Unit 1000

Pump Unit (UK) Order code: Med/P1000

Power: 230v +/-10% - 50Hz 40VA

Pump Unit (USA) Order code: Med/P1001

Power: 115v +/-10% - 60Hz 40VA

Pump Unit (EURO) Order code: Med/P1002

Power: 230v +/-10% 50Hz 40VA

Pump Unit Dimensions:

Height: 23.5 cm (9.25")

Width: 30 cm (12")

Depth: 16.6 cm (6.5")

Weight: 3.15 Kg (6.9 lbs)

Warranty: 2 Years on Pump, 1 Year on mattress.

Mattress Overlay 1000

Mattress only Order Code: Mat/1000ovm

Dimensions: 209 x 86 cm (81.5"x 34")

Air Cell Height Inflated: 12.7 cm (5.1")

Air Cells: 17 replaceable cells

Air Cell Material: PU (Polyurethane)

Cover Material: Nylon Coated Polyurethane

Patient Weight Limit: 140Kg (308lbs) (22 Stones)



Conforms to all relevant directives

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.



Medisphere

Global Wound Care Products

Medisphere Healthcare Ltd.
PO Box 346
High Wycombe.
Buckinghamshire.
HP13 7GD
United Kingdom.

e: inquires@medispherehealthcare.com

eurosales@medispherehealthcare.com

statesales@medispherehealthcare.com

t: +44 (0) 1494 522584

f: +44 (0) 1494 443468

w: www.medispherehealthcare.com

