


Medisphere introduces a new line in Wound Care Support Surfaces - 2000

The Tradewind 2000 mattress replacement system is the value for money, general-purpose dynamic support surface, providing effective pressure relief with durability and optimum Patient comfort.

The Tradewind 2000 is a true alternating support surface providing effective pressure relief by gently inflating and deflating the individual cells alternately over a ten minute cycle.

Patients Assessed at Medium to High Risk



www.medispherehealthcare.com

Technical Specifications
Tradewind 2000 Mattress Replacement System
 Complete System Order Code: Med/2000/MR

Pump Unit 2000

Pump Unit (UK) Order code: Med/P2000
 Power: 230v +/-10% - 50Hz 40VA

Pump Unit (USA) Order code: Med/P2001
 Power: 115v +/-10% - 60Hz 40VA

Pump Unit (EURO) Order code: Med/P2002
 Power: 230v +/-10% - 50Hz 40VA

Pump Unit Dimensions:
 Height: 23.5 cm (9.25")
 Width: 30 cm (12")
 Depth: 16.6 cm (6.5")
 Weight: 3.15 Kg (6.9 lbs)


Warranty: 2 Years on Pump, 1 Year on mattress.

Mattress Replacement 2000

Mattress only Order code: Mat/2000/MR
 Dimensions: 209 x 86 cm (81.5"x 34")
 Air Cell Height Inflated: 12.7 cm (5.1")
 Air Cells: 17 Replaceable cells
 Mattress Height Inflated: 20.3 cm (8")
 Foam Base: 3 inch
 Air Cell Material: PU (Polyurethane)
 Cover Material : Nylon Coated Polyurethane
 Patient Weight Limit: 160Kg (352lbs) (25 Stones)

CE
 Conforms to all relevant directives.

As our policy is one of continual improvement, we reserve the right to modify designs without prior notice.



Information on other Wound Care Support Surfaces from Medisphere Healthcare are available on request.

View our website for more information.

www.medispherehealthcare.com

Medisphere
 Global Wound Care Products

Medisphere Healthcare Ltd.
 PO Box 346
 High Wycombe.
 Buckinghamshire.
 HP13 7GD
 United Kingdom.

e: inquires@medispherehealthcare.com
eurosales@medispherehealthcare.com
statesales@medispherehealthcare.com
 t: +44 (0) 1494 522584
 f: +44 (0) 1494 443468
 w: www.medispherehealthcare.com