



www.medispherehealthcare.com

Medisphere introduces a new line in Wound Care Support Surfaces - 3000

Technical Specifications

Tradewind 3000 Mattress Replacement System

Complete System Order Code: Med/3000/MR

Pump Unit 3000

Pump Unit (UK) Order code: Med/P3000

Power: 230v +/-10% - 50Hz 40VA

Pump Unit (USA) Order code: Med/P3001

Power: 115v +/-10% - 60Hz 40VA

Pump Unit (EURO) Order code: Med/P3002

Power: 230v +/-10% - 50Hz 40VA

Pump Unit Dimensions:

Height: 23.5 cm (9.25")

Width: 30 cm (12")

Depth: 16.6 cm (6.5")

Weight: 3.15 Kg (6.9 lbs)

Warranty: 2 Years on Pump,1 Year on mattress.

Mattress Replacement 3000

Mattress only Order Code: Mat/3000/MR

Dimensions: 208.5 x 87.6 cm (82" x 34.5")

Air Cell Height Inflated: 20 cm (8")

Air Cell Material: PU (Polyurethane)

Air Cells: 20 Replaceable cells


Cover Material: Nylon Coated Polyurethane

Patient Weight Limit: 160 Kg (352lbs) (25 Stones)



Conforms to all relevant directives.

As our policy is one of continuous improvement,we reserve the right to modify designs without prior notice.



Medisphere

Global Wound Care Products

Medisphere Healthcare Ltd.
PO Box 346
High Wycombe.
Buckinghamshire.
HP13 7GD
United Kingdom.

e: inquires@medispherehealthcare.com

eurosales@medispherehealthcare.com

statesales@medispherehealthcare.com

t: +44 (0) 1494 522584

f: +44 (0) 1494 443468

w: www.medispherehealthcare.com