




Medisphere introduces a new line in Wound Care Support Surfaces - 3500

Tradewind 3500
 has been designed as a economy combination alternating low air loss therapy system with user friendly features, simplified everyday patient management. For the prevention and treatment of pressure sores in the medium to high risk patient.
Patients Assessed at Medium to High Risk



www.medispherehealthcare.com

Technical Specifications

*Tradewind 3500 Economy Combination
 Alternating with Low Air Loss Mattress Replacement System*

Complete System Order Code: Med/3500/lal MR

Pump Unit 3500
Pump Unit (UK) Order Code: Med/P3500
Power: 230v +/-10% - 50Hz 40VA
Pump Unit (USA) Order Code: Med/P3501
Power: 115v +/-10% - 60Hz 40VA
Pump Unit (EURO) Order Code: Med/P3502
Power: 230v +/-10% - 50Hz 40VA

Pump Unit Dimensions:
 Height: 23.5 cm (9.25")
 Width: 30 cm (12")
 Depth: 16.6 cm (6.5")
 Weight: 3.15 Kg (6.9lbs)



Warranty: 2 years on pump, 1 year on mattress.

Mattress Replacement 3500

Mattress Only Order Code: Mat/3500/lal MR
Dimensions: 208 x 87.6 cm (82"x 34.5")
Air Cell Height Inflated: 20 cm (8")
Air Cell Material: PU (Polyurethane)
Air Cells: 20 Replaceable Cells
Cover Material: Polyurethane Coated Nylon
Patient Weight Limit: 203Kg (446lbs) (32 Stones)



Conforms to all relevant directives.
 As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.



Information on other Wound Care Support Surfaces from Medisphere Healthcare are available on request.
 View our website for more information.
 www.medispherehealthcare.com

Medisphere
 Global Wound Care Products

Medisphere Healthcare Ltd.
 PO Box 346
 High Wycombe
 Buckinghamshire
 HP13 7GD
 United Kingdom

e: inquires@medispherehealthcare.com
eurosales@medispherehealthcare.com
statesales@medispherehealthcare.com
 t: +44 (0) 1494 522584
 f: +44 (0) 1494 443468
 w: www.medispherehealthcare.com