



Medisphere introduces a new line in Wound Care Support Surfaces - 3500

The Tradewind 3500 has been designed as an economy combination alternating low air loss therapy system with user friendly features, simplified everyday Patient management. For the prevention and treatment of pressure sores in the Medium to High Risk Patient.

Patients Assessed at Medium to High Risk

www.medispherehealthcare.com

Technical Specifications

*Tradewind 3500 Economy Combination
Alternating with low air loss Mattress Replacement System
Complete System Order Code: Med/3500/lal MR*

Pump Unit 3500

*Pump Unit (UK) Order code: Med/P3500
Power: 230v +/-10% - 50Hz 40VA*

*Pump Unit (USA) Order code: Med/P3501
Power: 115v +/-10% - 60Hz 40VA*

*Pump Unit (EURO) Order code: Med/P3502
Power: 230v +/-10% - 50Hz 40VA*

*Pump Unit Dimensions:
Height: 23.5 cm (9.25")
Width: 30 cm (12")
Depth: 16.6 cm (6.5")
Weight: 3.15 Kg (6.9 lbs)*

Warranty: 2 Years on Pump, 1 Year on mattress.


Mattress Replacement 3500

*Mattress Only Order Code: Mat/3500/lal MR
Dimensions: 208 x 87.6 cm (82"x 34.5")
Air Cell Height Inflated: 20 cm (8")
Air Cell Material: PU (Polyurethane)
Air Cells: 20 Replaceable cells
Cover Material: Polyurethane coated Nylon
Patient Weight Limit: 203Kg (446lbs) (32 Stones)*



Conforms to all relevant directives.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.



Medisphere
Global Wound Care Products

Medisphere Healthcare Ltd.
PO Box 346
High Wycombe.
Buckinghamshire.
HP13 7GD
United Kingdom.

e: inquires@medispherehealthcare.com
eurosales@medispherehealthcare.com
statesales@medispherehealthcare.com
t: +44 (0) 1494 522584
f: +44 (0) 1494 443468
w: www.medispherehealthcare.com