



Medisphere introduces a new line in Wound Care Support Surfaces - 5000

**The Tradewind 5000  
Combination Therapy Mattress  
Replacement System**

Provides Advanced Alternating Pressure  
Relief and Low Air Loss  
Two best therapies combined in one  
Deep Air Cell Cushion system.

*Patients Assessed at High to Very High Risk*

[www.medispherehealthcare.com](http://www.medispherehealthcare.com)

**Technical Specifications**

*Tradewind 5000 Alternating Low Air loss  
Deep Cell Mattress Replacement System*

*Complete System Order Code: Med/5000/lal MR*

**Pump Unit 5000**

*Pump Unit (UK) Order code: Med/P5000*

*Power: 230v +/-10% - 50Hz 80VA*

*Pump Unit (USA) Order code: Med/P5001*

*Power: 115v +/-10% - 60Hz 80VA*

*Pump Unit (EURO) Order code: Med/P5002*

*Power: 230v +/-10% 50Hz 80VA*

**Pump Unit Dimensions:**

*Height: 23.5 cm ( 9.25" )*

*Width: 30 cm (12" )*

*Depth: 16.6 cm ( 6.5" )*

*Weight: 3.15 Kg ( 6.9 lbs)*

*Remote Unit Order code: Remote/5000*

*Warranty: 2 Years on Pump, 1 Year on mattress.*

**Mattress Replacement 5000**

*Mattress only Order Code: Mat/lal 5000 MR*

*Dimensions: 208.5 x 87.6 cm (82"x 34.5")*

*Air Cell Height Inflated: 20 cm (8")*

*Air cell material: PU (Polyurethane)*

*Air cells: 20 replaceable cells*

*Cover Material: Nylon Coated Polyurethane*

*Patient Weight Limit: 286 Kg (308lbs) (45 Stones)*

*Bariatric Mattress Order Code: Bariatric/5000*



Conforms to all relevant directives

As our policy is one of continuous improvement,  
we reserve the right to modify designs without  
prior notice.




Medisphere 1000



Medisphere 2000



Medisphere 3000




Medisphere 600




Medisphere 3500



Medisphere 4000




Medisphere 7000



Medisphere 6000



Medisphere 5000



Information on other  
Wound Care Support Surfaces  
from  
Medisphere Healthcare  
are available on request.  
  
View our website  
for more information.  
[www.medispherehealthcare.com](http://www.medispherehealthcare.com)



**Medisphere**

*Global Wound Care Products*

Medisphere Healthcare Ltd.  
PO Box 346  
High Wycombe.  
Buckinghamshire.  
HP13 7GD  
United Kingdom.

e: [inquires@medispherehealthcare.com](mailto:inquires@medispherehealthcare.com)

[eurosales@medispherehealthcare.com](mailto:eurosales@medispherehealthcare.com)

[statesales@medispherehealthcare.com](mailto:statesales@medispherehealthcare.com)

t: +44 (0) 1494 522584

f: +44 (0) 1494 443468

w: [www.medispherehealthcare.com](http://www.medispherehealthcare.com)