


Medisphere introduces a new line in Wound Care Support Surfaces - 6000



**The Tradewind 6000**  
 Our most Advanced Combination Therapy System that incorporates new developments reflecting your needs within the ever changing acute and community care environments.

Alternating Pressure Relief / Low Air Loss therapy mattress replacement system for those Patients assessed at High to Very High Risk of developing Pressure Sores or with Superficial Established pressure sores.

*Patients Assessed at High to Very High Risk*



www.medispherehealthcare.com

**Technical Specifications**  
 Tradewind 6000 Advanced Dynamic with Low air Loss Mattress Replacement System  
 Complete System Order Code: Med/6000/ADL MR


**Pump Unit 6000**  
 Pump Unit (UK) Order code: Med/P6000  
 Power: 230v +/-10% - 50Hz 80VA  
 Pump Unit (USA) Order code: Med/P6001  
 Power: 115v +/-10% - 60Hz 80VA  
 Pump Unit (EURO) Order code: Med/P6002  
 Power: 230v +/-10% - 50Hz 80VA

**Pump Unit Dimensions:**  
 Height: 23.5 cm ( 9.25" )  
 Width: 30 cm (12")  
 Depth: 16.6 cm ( 6.5" )  
 Weight: 3.15 Kg ( 6.9 lbs)

Remote Unit Order Code: Remote/6000  
 Warranty: 2 Years on Pump, 1 Year on mattress.

**Mattress Replacement 6000**  
 Mattress only Order Code: Mat/ADL 6000/MR  
 Dimensions: 208.5 x 87.6 cm (82" x 34.5")  
 Air Cell Height Inflated: 20 cm (8")  
 Air Cell Material: PU (Polyurethane)  
 Cover Material: Nylon Coated Polyurethane  
 Patient Weight Limit: 286Kg (630lbs) (45 Stones)  
 Bariatric Mattress Order Code: mat-Bariat/6000

**CE**  
 Conforms to all relevant directives.  
 As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.



**Medisphere**  
 Global Wound Care Products

Medisphere Healthcare Ltd.  
 PO Box 346  
 High Wycombe.  
 Buckinghamshire.  
 HP13 7GD  
 United Kingdom.

e: [inquires@medispherehealthcare.com](mailto:inquires@medispherehealthcare.com)  
[eurosales@medispherehealthcare.com](mailto:eurosales@medispherehealthcare.com)  
[statesales@medispherehealthcare.com](mailto:statesales@medispherehealthcare.com)  
 t: +44 (0) 1494 522584  
 f: +44 (0) 1494 443468  
 w: [www.medispherehealthcare.com](http://www.medispherehealthcare.com)