



Medisphere introduces a new line in Wound Care Support Surfaces - 7000

Technical Specifications

Cyclone 7000 Three Zone Deep Cell Low Air Loss Mattress Replacement System

Complete System Order Code: Med/7000/z-lal MR

Pump Unit 7000

Pump Unit (UK) Order code: Med/P7000

Power: 230v +/-10% - 50Hz 405VA

Pump Unit (USA) Order code: Med/P7001

Power: 115v +/-10% - 60Hz 405VA

Pump Unit (EURO) Order code: Med/P7002

Power: 230v +/-10% - 50Hz 405VA

Pump Unit Dimensions:

Height: 23.5 cm (9.25")

Width: 30 cm (12")

Depth: 16.6 cm (6.5")

Weight: 3.15 Kg (6.9 lbs)

Warranty: 2 Years on Pump, 1 Year on mattress.

Mattress Replacement 7000

Mattress only Order Code: Mat/7000/MR

Dimensions: 208.5 x 87.6 cm (82" x 34.5")

Air Cell Height Inflated: 20 cm (8")

Air Cell Material: PU (Polyurethane)

Cover Material: Nylon Coated Polyurethane

Patient Weight Limit: 203Kg (446lbs) (32 Stones)

Bariatric Mattress Order Code: mat-Bariat/7000



Conforms to all relevant directives.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.




Information on other Wound Care Support Surfaces from Medisphere Healthcare are available on request. View our website for more information. www.medispherehealthcare.com

Medisphere
Global Wound Care Products

Medisphere Healthcare Ltd.
PO Box 346
High Wycombe.
Buckinghamshire.
HP13 7GD
United Kingdom.

e: inquires@medispherehealthcare.com
eurosales@medispherehealthcare.com
statesales@medispherehealthcare.com
t: +44 (0) 1494 522584
f: +44 (0) 1494 443468
w: www.medispherehealthcare.com

The Cyclone 7000 Low Air Loss Mattress Replacement System provides excellent pressure reduction for Patients assessed up to Very High risk of developing pressure sores, and including group 4. Users unable to tolerate a moving pressure relieving support surface may benefit from Low Air Loss therapy. "Simplicity, Reliability and Value in One" Patients Assessed at High to Very High Risk



www.medispherehealthcare.com